


Guardian Call

The GuardianCall system offers real-time notifications, medication reminders, and SOS alerts to support caretakers and relatives. It includes a main device, an SOS button, and a mobile app, ensuring your loved ones are always connected and cared for.

This device Has 3 main components

- Main Device
- SOS Button
- Mobile App



- Main Device
- SOS Button




- Mobile App




- Pen for Easy touch

SOS Button



- Charging Port: Charges using a standard USB-C charger.
- Pressing Button: sends an emergency alert to all the caretaker's smartphones.
- Strap is put around the wrist
- One battery Cycle last upto 20 days

Main Device



For Charging main device

For charging the SOS button

1. Initial Setup:

- On first boot, select a Wi-Fi network and enter the password using the provided PIN. The password won't be needed again for the same network.

2. Wi-Fi Connection:

- The device auto-connects if the network is available. If Wi-Fi is off, a pop-up will prompt you to switch networks or retry. Pressing space registers a new Wi-Fi, which will auto-connect once available.

3. No Wi-Fi Connection:


- The device beeps continuously if no Wi-Fi is available. Connect to Wi-Fi to use the mobile app.

4. Charging:

- Charge the device via the USB-C port. Use the USB-A port to charge the button.


Mobile App

1. Login and Register



- Once the app is open, you will be directed to the login page
- Tap on the field and enter your registered username or email address.

If you haven't registered yet, click register



- **Step 2: Enter Device Credentials**
Apart from the usual email and password fill in the following fields



Device ID:

- Enter the unique Device ID provided to you at the time of purchase.

Secret Code:



- Enter the Secret Code given to you with your device. This code is case-sensitive, so make sure to enter it correctly.

2. Buttons




- To initially display a button in the app, the user must have pressed the button at least once
- the user can rename the buttons that are displayed on the app as they see fit. the names may differ from user to user.

3. Setting Alarms




- To set an alarm, select the time, choose the day/days of the week, and enter a custom message for the alarm.
- Tap "save" to set the alarm



- Tap the bell icon to set a new alarm
- Tap to view the Alarm message or to remove existing alarms

3. Removing registered caretaker accounts



- Select no longer needed registered caretaker accounts
- Tap "remove" to remove the selected accounts



- When a patient is having an emergency a notification will be displayed along with a ringing tone
- when the notification is clicked this screen would be displayed
- Tap okay to stop the app from ringing,